

SFWRENG 3K04

Part 2 Documentation - DCM Design

Group #12

Matt Mione, Riley Mione, Derek Paylor, Johnathan Spinelli, Laura Yang

Table of Contents

1.1 Module: dcm_main.py	2
1.1.1 Purpose of Module	2
1.1.2 Black Box Behaviour of Module	2
1.1.3 Secret of Module	3
1.2 Module: Interface.py	3
1.2.1 Purpose of Module	3
1.2.2 Global Variables in Module	3
1.2.3 Private Functions in Module	3
1.2.4 Internal Behaviour of Each Function	4
1.2.5 Module's Secret	5
1.3 Module: MiscFunctions.py	5
1.3.1 Purpose of Module	5
1.3.2 Global Variables in Module	5
1.3.3 Constants	5
1.3.4 Public Functions Provided by Module	6
1.3.5 Private Functions in Module	6
1.3.6 Internal Behaviour of Each Function	6
1.3.7 Module's Secret	7
2. Requirements Likely to Change	8
3. Design Decisions Likely to Change	8
4. Testing DCM	9
4.1 User Creation Procedure:	9
4.2 Login Procedure:	10
4.3 Program Interface:	12
4.4 Rate-Adaptive Modes:	37
4.5 Improper Transmission/Other Issues	38

1.1 Module: dcm_main.py

1.1.1 Purpose of Module

This module is designed to be the main component of the program. Akin to an executable, this module is used to “launch” the program. With minimal code, this module interacts with the other components of the program.

1.1.2 Black Box Behaviour of Module

Upon launching the program, a small Welcome screen is presented to the user. There are two buttons, one says “Login” and the other says “Create User”. By clicking the Login button, the screen is changed to a screen with a message saying “Enter your credentials”, two entry boxes, a “Submit” button and a “Return to start page” button. Clicking the Return button takes the user back to the Welcome screen. If the user clicks the Submit button with no credentials or invalid credentials typed into the entry boxes, a pop-up window informs the user that they entered invalid credentials. If the user enters valid credentials, a new window appears (see next paragraph). If the user had clicked on the Create User button from the main menu, they would be taken to a similar looking screen as the Login screen, but this one says “Create a new user” at the top. Clicking the Return button takes the user back to the Welcome screen. If the user clicks the Submit button without entering a username, a pop-up telling the user to input a username. If the user clicks the Submit button with an existing username in the entry box, a pop-up tells the user to login with that existing account, or enter a different username. If the user clicks Submit without a password entered, a pop-up tells the user to enter a password. If the user clicks the Submit button with a unique username and a valid password, a new window appears (see next paragraph) and a new entry in the Database.JSON file appears.

The new window contains 11 tabs at the top, listing the different modes for the pacemaker to be programmed in. Each tab contains a varying number of number sliders with different labels beside each. At the bottom of the screen, the user’s username is displayed. There is a Submit button, as well as three buttons that say “Atrium Egram”, “Ventricle Egram” and “Dual Egram”. When the Atrium Egram button is clicked, a line graph with a blue line pops up. When the Ventricle Egram button is clicked, a line graph with a red line pops up. When the Dual Egram button is clicked, a line graph containing both the red line and blue line appears. The user is able to interact with the sliders, which changes the number above the slider bar. Upon clicking the submit button, a pop-up tells the user that the Pacemaker values were updated successfully. If the “Lower Rate Limit” value is greater than the “Upper Rate Limit” value, the pop-up also says that the Lower Rate Limit was fixed to 50ppm. Another pop-up tells the user that the database was successfully updated with their credentials. The file called “SerialComm.JSON” is updated with new values. The username of the current user is written in the SerialComm file, with the rest of the programmable parameters after it. If the parameter name appeared on the tab that the user clicked submit from, that parameter has the same value in the SerialComm file, but even if that parameter did not appear on that tab, the parameter still appears in the file with a value. If the pacemaker is connected to the computer upon clicking the submit button, the pacemaker will behave according to the programmed parameters.

1.1.3 Secret of Module

The secret of *dcm_main.py* is that it utilizes the other two modules of the program, *Interface.py* and *MiscFunctions.py*, to do the brunt of the work for it. The workflow is split up for neatness and separation of concerns for the program. *Interface* handles the GUI portion and most of the visual aspects of the DCM. It calls upon *MiscFunctions* for many of its processes, such as dumping and encrypting data, updating the program's user database, and the pacemaker parameters for each user, as well as for handling the serial communication between the microcontroller and the DCM host machine. Essentially, *dcm_main.py* is the main executable file, so-to-speak for the entire application.

1.2 Module: Interface.py

1.2.1 Purpose of Module

The interface module is designed to handle everything to do with user interaction. This module contains classes and functions that manage the menus and navigation through the program.

1.2.2 Global Variables in Module

This module contains three global variables.

UserID: this global variable contains the current user's username, for use in several places in the module.

AEgram: this variable was one of two to be used for Egram processing, but ended up not being used, due to the requirements change that Egram will not be mandatory.

VEgram: the other of two variables to be used for Egram processing, but not currently in use for reasons stated above.

1.2.3 Private Functions in Module

def _switch_frame(frame_class):

Takes a Tkinter frame object and switches to that part of the interface.

def _create_user(username, password, frame_class)

This function takes the desired username and password input into the two text boxes in the GUI and does some operations. It checks it against the existing database, which must be limited to 10 users, to see if it exists, and if not, will input it into the database. Also takes a frame_class object.

def _login_test(username, password, frame_class)

Function takes text-variable version of username and password, as well as frame class.

def _egramSwitch(value)

Egram has not been implemented fully yet.

def _animate(i)

Egram has not been implemented fully yet.

1.2.4 Internal Behaviour of Each Function**def _switch_frame(frame_class):**

This is a function contained in the DCM class. It takes the current Tkinter frame, and then replaces it with the inputted frame_class object. It does this through using the tkinter frame.destroy() function and then setting the window's frame equal to the desired frame.

def _create_user(username, password, frame_class)

This function encrypts the given username and password using the encrypt() function. After encryption, the encrypted values are compared against the existing database of users. In order for user creation to be successful, the database must have less than 10 existing users, the username must not already exist, and at least one character has to be typed for both username and password. If these conditions are met, the database is updated with the user's encrypted credentials, and the current frame is switched. If one or more of the conditions are not met, a message box pops up with the first condition that was not met.

def _login_test(username, password, frame_class)

Function encrypts the given username and password combination. If the encrypted username and password do not match the database, the combination is rejected and the user is informed of the condition that was not met. If the username does not exist, the user is also informed.

def _egramSwitch(value)

This function controls which set of ECG data would be displayed for the user. Depending on the value that is passed to the function, the global variables AEgram and VEgram are altered to define which values should be displayed on the plot. Calls animation function.

def _animate(i)

Reads values from text files and plots values on the plot. Depending on the values of AEgram and VEgram, it selects which text file to read from.

1.2.5 Module's Secret

To prevent the user from inputting invalid pacemaker parameters, the interface uses sliders to obtain values. The sliders are able to limit the input values between a maximum and minimum acceptable value

with appropriate increments. To reduce the number of functions in the module, the same function is used to collect values from each parameter tab. When a programmable parameter isn't used in a specific mode, the program passes that parameter's "nominal value", which is a value within its functioning range.

1.3 Module: MiscFunctions.py

1.3.1 Purpose of Module

This module controls the data interactions and other miscellaneous functions. This involves the databases where the data is stored and communication of data between devices. It contains many public functions used by other modules as well.

1.3.2 Global Variables in Module

database: dictionary used to store all username and password combinations of existing users.

pacemaker_values: dictionary used to store pacemaker parameters for each user.

board: a PySerial object used to communicate over USB with the microcontroller.

1.3.3 Constants

SHIFT: constant value used for encryption

baud_rate: constant value used for serial communication

UPLOAD_LOCATION: constant string variable that is declared globally, denotes where encrypted user database will be dumped to on local machine.

DUMP_LOCATION: constant string variable that is declared globally, denotes where the pacemaker's inputted values will be dumped in the local machine.

```
''' DATABASE '''
# Checks if the json file exists
if os.path.exists(DUMP_LOCATION):
    # if it exists, load in that database as the current database. Now it has memory
    with open(DUMP_LOCATION) as f:
        database = json.load(f)
```

SHIFT: constant integer used in encryption and decryption processes.

1.3.4 Public Functions Provided by Module

def update_info(mode, low, up, AAmp, VAmp, APW, VPW, ASense, VSense, ARP, VRP, MaxSense, PVARP, FAVD, ReTime, RecTime, RespFact, AThresh, user)

Neatly updates dictionary with pacemaker parameters as per requirements in documentation and communicates the parameters to the pacemaker.

class IO (*args, **kwargs)

Simple class to do with file i/o and encryption/decryption. Consists of three class methods.

def decrypt(some_phrase)

This function takes the given phrase and returns the decrypted message.

def dump(path, data_dict)

This function writes the current state of the user database to a json file. Having a local copy is very important, since this copy can be written to memory and accessed whether or not the DCM python script is running or not, and when the script reboots it remembers who has been registered.

def encrypt(some_phrase)

This function takes the given phrase and returns the encrypted message.

def communicate_parameters(mode, low, up, AAmp, VAmp, APW, VPW, ASense, VSense, ARP, VRP, MaxSense, FAVD, ReTime, RecTime, RespFact, AThresh):

Takes input values and communicates them to the board.

1.3.5 Private Functions in Module

def _to_bytes(mode, low, up, Aamp, Vamp, Apw, Vpw, Asense, Vsense, ARP, VRP, MSR, FAVD, RE, REC, RES, AT):

Takes input variables and packs them into a byte array, including parity byte at the end.

1.3.6 Internal Behaviour of Each Function

def update_info(mode, low, up, AAmp, VAmp, APW, VPW, ASense, VSense, ARP, VRP, MaxSense, PVARP, FAVD, ReTime, RecTime, RespFact, AThresh, user)

This function checks the parameters *low* and *up*, the lower and upper pacing thresholds respectively, and makes sure that *low* is not greater than *up*, and if it is, *low* is set equal to 50 to circumvent that from happening. If the *mode* involves a maximum sensor rate, the function

checks that *low* is less than the maximum sensor rate. If it isn't, the function returns and displays a pop-up that informs the user of this issue. If there is no maximum sensor rate issue, the program writes the pacemaker values for that username to the JSON file *SerialComm.json*. The function uses the *communicate_parameters* function to send the parameters to the pacemaker.

class IO (*args, **kwargs)

def decrypt(some_phrase)

This function uses the constant SHIFT to negatively shift the ASCII value of each letter, which effectively decrypts the phrase.

def dump(path, data_dict)

The contents in *data_dict* are written into the path and formats it in a way that is easy to read, both for the user and the program.

def encrypt(some_phrase)

This function uses the constant SHIFT to positively shift the ASCII value of each letter, which effectively encrypts the phrase.

def communicate_parameters(mode, low, up, AAmp, VAmp, APW, VPW, ASense, VSense, ARP, VRP, MaxSense, FAVD, ReTime, RecTime, RespFact, AThresh):

Receives input parameters and converts them to a format that can be used to communicate to the board. Sends the parameters to the board and confirms that the board accepted the values. Compares parameters on the board with parameters from DCM to make sure there are no inconsistencies.

def _to_bytes(mode, low, up, Aamp, Vamp, Apw, Vpw, Asense, Vsense, ARP, VRP, MSR, FAVD, RE, REC, RES, AT):

Converts each programmable parameter into the data type accepted by the pacemaker and forms them into a byte array with a parity byte at the end.

1.3.7 Module's Secret

This module uses a very basic encryption method to keep user credentials safe. This is used to make sure that if someone were to read the database file, they wouldn't know every user's credentials. By simply changing the ASCII value of each character, the credentials are encrypted and decrypted.

To communicate with the pacemaker, the DCM uses serial communication. The PySerial object, *board*, needs to be configured for the host machine. The *port* parameter needs to be changed to the USB port that

the microcontroller is plugged into. Serial communication allows values to be transmitted very quickly to the pacemaker. To ensure that the values were being passed in the correct order and were passed correctly, the pacemaker sends the received values back to the DCM to confirm.

2. Requirements Likely to Change

The initial requirements for the DCM specified a functional ECG component. The requirements have changed and no longer require an ECG. Development for displaying the ECG was started and included, but it is not functional. For future models of the DCM, the ECG would likely be implemented and functioning, or removed altogether. The initial requirements also required two-way communication between the board and the DCM, but communication from the board to the DCM is no longer required. Two-way communication is already implemented, but not to the point where it would work for the ECG. Future models will have communication implemented for more advanced communication between the board and the DCM.

3. Design Decisions Likely to Change

The program is split into five different files - 3 python files and 2 JSON files. The python files are separated based on functionality, but they could be separated further based on their content. The program functions properly with the current configuration, but to make it easier to adjust or implement new functions, the files can be made more specific.

The GUI consists of multiple windows with varying sizes. The largest window takes up the majority of the screen and would not be ideal for a device with a small screen size. A future implementation would likely have everything be displayed on a single window. For accessibility purposes, the window could have a scrollbar that allows users to navigate through the pages with ease.


Modularity was a large factor in the design of the program. Separating the program into modules made development easy, as team members were able to work on their own module without any conflict. Modularity allowed for separation of concerns. Components that deal with communicating with the pacemaker and data files are encapsulated within one module, while other functions dealing with the interface are encapsulated in another. This created high cohesion and low coupling in the program.


In the case where an external element is missing (ie. file doesn't exist, file name change), the program allows for an easy correction. File names and databases are declared and visible near the top of each module. Data access or delivery is broken down into small functions, that make changing the code easy. For example, the ECG functionality is not fully implemented, but since the process is composed of smaller processes, the code can be altered easily to change how the values are stored and accessed.

4. Testing DCM

4.1 User Creation Procedure:


A few errors can occur in the user creation procedure. These being if the database is full, if the username exists, and if there is no password or username entered and submit is clicked.






4.2 Login Procedure:


Here we entered a username that didn't exist, expecting it to show a User Does Not Exist window, and it had the expected behaviour.




Here we entered the incorrect password, expecting it to show that the password was incorrect. It behaved as expected.



This case involves having either the username or password be blank. The program performed as expected and sent out the corresponding window.



Here we are testing the case that the user exists and that the password is correct. The login should go smoothly.



4.3 Program Interface:

Proper Saves and Transmissions

We will need to test the submitting of parameters for each of the modes. Here we submit good values for all of them. The values transmit and save properly, as expected of them.

AOO

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDR

Lower Rate Limit (ppm) 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 90 105 120 135 150 165

Atrial Amplitude(V) 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Atrial Pulse Width (ms) 1 2 3 4 5 6 7 8 9 10

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDR

Lower Rate Limit (ppm) 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 90 105 120 135 150 165

Atrial Amplitude(V) 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Atrial Pulse Width (ms) 1 2 3 4 5 6 7 8 9 10

Submit


User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

All parameters were transmitted properly

OK

**VOO**

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 90 105 120 135 150 165

Ventricular Amplitude (V) 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Ventricular Pulse Width (ms) 1 2 3 4 5 6 7 8 9 10

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

All parameters were transmitted properly

OK

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 90 105 120 135 150 165

Ventricular Amplitude (V) 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Ventricular Pulse Width (ms) 1 2 3 4 5 6 7 8 9 10

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Message

Pacemaker Values Updated Successfully

OK

AAI

Pacemaker Interface

AOO VOO **AAI** VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 90 105 120 135 150 165

Atrial Amplitude (V) 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Atrial Pulse Width (ms) 1 2 3 4 5 6 7 8 9 10

Atrial Sense Threshold (V) 0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2

ARP (ms) 150 170 190 210 230 250 270 290 310 330 350 370 390 410 430 450 470 490

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

All parameters were transmitted properly

OK

AOOR

Pacemaker Interface — □ ×

AOO VOO AAI VVI DOO **AOOR** VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) <input type="text" value="50"/> <input type="text" value="50"/> 50 65 80 95 110 125 140 155 170	Atrial Pulse Width (ms) <input type="text" value="1"/> <input type="text" value="1"/> 1 2 3 4 5 6 7 8 9 10
Upper Rate Limit (ppm) <input type="text" value="75"/> <input type="text" value="75"/> 75 90 105 120 135 150 165	Reaction Time (s) <input type="text" value="1"/> <input type="text" value="1"/> 1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Max Sensor Rate (ppm) <input type="text" value="50"/> <input type="text" value="50"/> 50 60 70 80 90 100 110 120 130 140 150 160 170	Recovery Time(s) <input type="text" value="5"/> <input type="text" value="5"/> 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V) <input type="text" value="0.5"/> <input type="text" value="0.5"/> 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0	Response Factor (Slow-Fast) <input type="text" value="1"/> <input type="text" value="1"/> 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High) <input type="text" value="1"/> <input type="text" value="1"/> 1 2 3 4	
<div> <input type="button" value="Submit"/> </div>	
User: zxc	<div> <input type="button" value="Atrium Egram"/> <input type="button" value="Ventricle Egram"/> <input type="button" value="Dual Egram"/> </div>

Pacemaker Interface

AOO VOO AAI VVI DOO **AOOR** VOOR AAIR VVIR DOOR DDRR


Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="65"/> <input type="text" value="80"/> <input type="text" value="95"/> <input type="text" value="110"/> <input type="text" value="125"/> <input type="text" value="140"/> <input type="text" value="155"/> <input type="text" value="170"/>	Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/> <input type="text" value="90"/> <input type="text" value="105"/> <input type="text" value="120"/> <input type="text" value="135"/> <input type="text" value="150"/> <input type="text" value="165"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="3"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="60"/> <input type="text" value="70"/> <input type="text" value="80"/> <input type="text" value="90"/> <input type="text" value="100"/> <input type="text" value="110"/> <input type="text" value="120"/> <input type="text" value="130"/> <input type="text" value="140"/> <input type="text" value="150"/> <input type="text" value="160"/> <input type="text" value="170"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/> <input type="text" value="11"/> <input type="text" value="12"/> <input type="text" value="13"/> <input type="text" value="14"/> <input type="text" value="15"/> <input type="text" value="16"/>
Activity Threshold (Low-High)		<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/>	

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

 All parameters were transmitted properly

OK

Pacemaker Interface

AOO VOO AAI VVI DOO **AOOR** VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm)	50	Atrial Pulse Width (ms)	1
	50 65 80 95 110 125 140 155 170		1 2 3 4 5 6 7 8 9 10
Upper Rate Limit (ppm)	75	Reaction Time (s)	1
	75 90 105 120 135 150 165		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Max Sensor Rate (ppm)	50	Recovery Time(s)	5
	50 60 70 80 90 100 110 120 130 140 150 160 170		5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V)	0.5	Response Factor (Slow-Fast)	1
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)		1	
		1 2 3 4	
Submit			
User: zxc	Atrium Egram Ventricle Egram Dual Egram		

Pacemaker Message

i Pacemaker Values Updated Successfully

OK

VOOR

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR **VOOR** AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="65"/> <input type="text" value="80"/> <input type="text" value="95"/> <input type="text" value="110"/> <input type="text" value="125"/> <input type="text" value="140"/> <input type="text" value="155"/> <input type="text" value="170"/>	Ventricular Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/> <input type="text" value="90"/> <input type="text" value="105"/> <input type="text" value="120"/> <input type="text" value="135"/> <input type="text" value="150"/> <input type="text" value="165"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="3"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="60"/> <input type="text" value="70"/> <input type="text" value="80"/> <input type="text" value="90"/> <input type="text" value="100"/> <input type="text" value="110"/> <input type="text" value="120"/> <input type="text" value="130"/> <input type="text" value="140"/> <input type="text" value="150"/> <input type="text" value="160"/> <input type="text" value="170"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Ventricular Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/> <input type="text" value="11"/> <input type="text" value="12"/> <input type="text" value="13"/> <input type="text" value="14"/> <input type="text" value="15"/> <input type="text" value="16"/>
Activity Threshold (Low-High)		<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/>	


User: zxc

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR **VOOR** AAIR VVIR DOOR DDRR

Lower Rate Limit (ppm)	50	Ventricular Pulse Width (ms)	1
	50 65 80 95 110 125 140 155 170		1 2 3 4 5 6 7 8 9 10
Upper Rate Limit (ppm)	75	Reaction Time (s)	1
	75 90 105 120 135 150 165		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Max Sensor Rate (ppm)	50	Recovery Time(s)	5
	50 60 70 80 90 100 110 120 130 140 150 160 170		5 7 9 11 13 15 17 19 21 23 25 27 29
Ventricular Amplitude (V)	0.5	Response Factor (Slow-Fast)	1
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)		1	
		1	2 3 4
<input type="button" value="Submit"/>			
User: zxc		<input type="button" value="Atrium Egram"/> <input type="button" value="Ventricle Egram"/> <input type="button" value="Dual Egram"/>	

Communication OK

 All parameters were transmitted properly

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) 50 50 65 80 95 110 125 140 155 170

Upper Rate Limit (ppm) 75 75 90 105 120 135 150 165

Max Sensor Rate (ppm) 50 50 60 70 80 90 100 110 120 130 140 150 160 170

Ventricular Amplitude (V) 0.5 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0

Ventricular Pulse Width (ms) 1 1 2 3 4 5 6 7 8 9 10

Reaction Time (s) 1 1 3 5 7 9 11 13 15 17 19 21 23 25 27 29

Recovery Time(s) 5 5 7 9 11 13 15 17 19 21 23 25 27 29

Response Factor (Slow-Fast) 1 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

Activity Threshold (Low-High) 1 1 2 3 4

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Message

Pacemaker Values Updated Successfully

OK

AAIR

Pacemaker Interface

☐ AOO
 ☐ VOO
 ☐ AAI
 ☐ VVI
 ☐ DOO
 ☐ AOOR
 ☐ VOOR
 ☒ AAIR
 ☐ VVIR
 ☐ DOOR
 ☐ DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	ARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	50 65 80 95 110 125 140 155 170		150 180 210 240 270 300 330 360 390 420 450 480
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/>	PVARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	75 90 105 120 135 150 165		150 180 210 240 270 300 330 360 390 420 450 480
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/>
	50 60 70 80 90 100 110 120 130 140 150 160 170		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/>
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Atrial Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Activity Threshold (Low-High)	<input type="text" value="1"/> <input type="text" value="1"/>
	0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2		1 2 3 4

Submit

User: zxc

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR **AAIR** VVIR DOOR DDDR


Lower Rate Limit (ppm)	50	ARP (ms)	150
<input type="text"/>	50 65 80 95 110 125 140 155 170	<input type="text"/>	150 180 210 240 270 300 330 360 390 420 450 480
Upper Rate Limit (ppm)	75	PVARP (ms)	150
<input type="text"/>	75 90 105 120 135 150 165	<input type="text"/>	150 180 210 240 270 300 330 360 390 420 450 480
Max Sensor Rate (ppm)	50	Reaction Time (s)	1
<input type="text"/>	50 60 70 80 90 100 110 120 130 140 150 160 170	<input type="text"/>	1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V)	0.5	Recovery Time(s)	5
<input type="text"/>	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0	<input type="text"/>	5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Pulse Width (ms)	1	Response Factor (Slow-Fast)	1
<input type="text"/>	1 2 3 4 5 6 7 8 9 10	<input type="text"/>	1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Atrial Sense Threshold (V)	0.5	Activity Threshold (Low-High)	1
<input type="text"/>	0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2	<input type="text"/>	1 2 3 4

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

 All parameters were transmitted properly

OK

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR **AAIR** VVIR DOOR DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	ARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	50 65 80 95 110 125 140 155 170		150 180 210 240 270 300 330 360 390 420 450 480
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/>	PVARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	75 90 105 120 135 150 165		150 180 210 240 270 300 330 360 390 420 450 480
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/>
	50 60 70 80 90 100 110 120 130 140 150 160 170		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/>
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Atrial Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Activity Threshold (Low-High)	<input type="text" value="1"/> <input type="text" value="1"/>
	0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2		1 2 3 4

Submit

User: zxc

Pacemaker Message

Pacemaker Values Updated Successfully

OK

gram Ventricle Egram Dual Egram

VVIR

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR **VVIR** DOOR DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="65"/> <input type="text" value="80"/> <input type="text" value="95"/> <input type="text" value="110"/> <input type="text" value="125"/> <input type="text" value="140"/> <input type="text" value="155"/> <input type="text" value="170"/>	Ventricular Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="0.8"/> <input type="text" value="1.1"/> <input type="text" value="1.4"/> <input type="text" value="1.7"/> <input type="text" value="2.0"/> <input type="text" value="2.3"/> <input type="text" value="2.6"/> <input type="text" value="2.9"/> <input type="text" value="3.2"/>
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/> <input type="text" value="90"/> <input type="text" value="105"/> <input type="text" value="120"/> <input type="text" value="135"/> <input type="text" value="150"/> <input type="text" value="165"/>	VRP (ms)	<input type="text" value="150"/> <input type="text" value="150"/> <input type="text" value="180"/> <input type="text" value="210"/> <input type="text" value="240"/> <input type="text" value="270"/> <input type="text" value="300"/> <input type="text" value="330"/> <input type="text" value="360"/> <input type="text" value="390"/> <input type="text" value="420"/> <input type="text" value="450"/> <input type="text" value="480"/>
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="60"/> <input type="text" value="70"/> <input type="text" value="80"/> <input type="text" value="90"/> <input type="text" value="100"/> <input type="text" value="110"/> <input type="text" value="120"/> <input type="text" value="130"/> <input type="text" value="140"/> <input type="text" value="150"/> <input type="text" value="160"/> <input type="text" value="170"/>	Reaction Time (s)	<input type="text" value="2"/> <input type="text" value="2"/> <input type="text" value="4"/> <input type="text" value="6"/> <input type="text" value="8"/> <input type="text" value="10"/> <input type="text" value="12"/> <input type="text" value="14"/> <input type="text" value="16"/> <input type="text" value="18"/> <input type="text" value="20"/> <input type="text" value="22"/> <input type="text" value="24"/> <input type="text" value="26"/> <input type="text" value="28"/> <input type="text" value="30"/>
Ventricular Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Ventricular Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/> <input type="text" value="11"/> <input type="text" value="12"/> <input type="text" value="13"/> <input type="text" value="14"/> <input type="text" value="15"/> <input type="text" value="16"/>
Activity Threshold (Low-High)		<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/>	

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR **VVIR** DOOR DDDR


Lower Rate Limit (ppm)	<input type="range" value="50"/>	Ventricular Sense Threshold (V)	<input type="range" value="0.5"/>
	50 65 80 95 110 125 140 155 170		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Upper Rate Limit (ppm)	<input type="range" value="75"/>	VRP (ms)	<input type="range" value="150"/>
	75 90 105 120 135 150 165		150 180 210 240 270 300 330 360 390 420 450 480
Max Sensor Rate (ppm)	<input type="range" value="50"/>	Reaction Time (s)	<input type="range" value="2"/>
	50 60 70 80 90 100 110 120 130 140 150 160 170		2 4 6 8 10 12 14 16 18 20 22 24 26 28 30
Ventricular Amplitude (V)	<input type="range" value="0.5"/>	Recovery Time(s)	<input type="range" value="5"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		5 7 9 11 13 15 17 19 21 23 25 27 29
Ventricular Pulse Width (ms)	<input type="range" value="1"/>	Response Factor (Slow-Fast)	<input type="range" value="1"/>
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)		<input type="range" value="1"/>	
		1 2 3 4	

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Communication OK

 All parameters were transmitted properly

OK

Pacemaker Interface

AOO | VOO | AAI | VVI | DOO | AOOR | VOOR | AAIR | **VVIR** | DOOR | DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/>	Ventricular Sense Threshold (V)	<input type="text" value="0.5"/>
	50 65 80 95 110 125 140 155 170		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Upper Rate Limit (ppm)	<input type="text" value="75"/>	VRP (ms)	<input type="text" value="150"/>
	75 90 105 120 135 150 165		150 180 210 240 270 300 330 360 390 420 450 480
Max Sensor Rate (ppm)	<input type="text" value="50"/>	Reaction Time (s)	<input type="text" value="2"/>
	50 60 70 80 90 100 110 120 130 140 150 160 170		2 4 6 8 10 12 14 16 18 20 22 24 26 28 30
Ventricular Amplitude (V)	<input type="text" value="0.5"/>	Recovery Time(s)	<input type="text" value="5"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		5 7 9 11 13 15 17 19 21 23 25 27 29
Ventricular Pulse Width (ms)	<input type="text" value="1"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/>
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)		<input type="text" value="1"/>	
		1 2 3 4	
<input type="button" value="Submit"/>			
User: zxc		<input type="button" value="Atrium Egram"/> <input type="button" value="Ventricle Egram"/> <input type="button" value="Dual Egram"/>	

Pacemaker Message

Pacemaker Values Updated Successfully

DOOR

Pacemaker Interface — □ ×

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR **DOOR** DDDR

Lower Rate Limit (ppm) <input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="65"/> <input type="text" value="80"/> <input type="text" value="95"/> <input type="text" value="110"/> <input type="text" value="125"/> <input type="text" value="140"/> <input type="text" value="155"/> <input type="text" value="170"/>	Atrial Pulse Width (ms) <input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>
Upper Rate Limit (ppm) <input type="text" value="75"/> <input type="text" value="75"/> <input type="text" value="90"/> <input type="text" value="105"/> <input type="text" value="120"/> <input type="text" value="135"/> <input type="text" value="150"/> <input type="text" value="165"/>	Ventricular Pulse Width (ms) <input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>
Max Sensor Rate (ppm) <input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="60"/> <input type="text" value="70"/> <input type="text" value="80"/> <input type="text" value="90"/> <input type="text" value="100"/> <input type="text" value="110"/> <input type="text" value="120"/> <input type="text" value="130"/> <input type="text" value="140"/> <input type="text" value="150"/> <input type="text" value="160"/> <input type="text" value="170"/>	Reaction Time (s) <input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="3"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Fixed AV Delay (ms) <input type="text" value="70"/> <input type="text" value="70"/> <input type="text" value="90"/> <input type="text" value="110"/> <input type="text" value="130"/> <input type="text" value="150"/> <input type="text" value="170"/> <input type="text" value="190"/> <input type="text" value="210"/> <input type="text" value="230"/> <input type="text" value="250"/> <input type="text" value="270"/> <input type="text" value="290"/>	Recovery Time(s) <input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Atrial Amplitude (V) <input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Response Factor (Slow-Fast) <input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/> <input type="text" value="11"/> <input type="text" value="12"/> <input type="text" value="13"/> <input type="text" value="14"/> <input type="text" value="15"/> <input type="text" value="16"/>
Ventricular Amplitude (V) <input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Activity Threshold (Low-High) <input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/>

User: zxc

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR


Lower Rate Limit (ppm)	50	Atrial Pulse Width (ms)	1
<input type="text"/>	50 65 80 95 110 125 140 155 170	<input type="text"/>	1 2 3 4 5 6 7 8 9 10
Upper Rate Limit (ppm)	75	Ventricular Pulse Width (ms)	1
<input type="text"/>	75 90 105 120 135 150 165	<input type="text"/>	1 2 3 4 5 6 7 8 9 10
Max Sensor Rate (ppm)	50	Reaction Time (s)	1
<input type="text"/>	50 60 70 80 90 100 110 120 130 140 150 160 170	<input type="text"/>	1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Fixed AV Delay (ms)	70	Recovery Time(s)	5
<input type="text"/>	70 90 110 130 150 170 190 210 230 250 270 290	<input type="text"/>	5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Amplitude (V)	0.5	Response Factor (Slow-Fast)	1
<input type="text"/>	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0	<input type="text"/>	1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Ventricular Amplitude (V)	0.5	Activity Threshold (Low-High)	1
<input type="text"/>	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0	<input type="text"/>	1 2 3 4

Submit

User: zxc

Atrium Egram Ventricule Egram Dual Egram

Communication OK

 All parameters were transmitted properly

OK

Pacemaker Interface


☐ AOO
 ☐ VOO
 ☐ AAI
 ☐ VVI
 ☐ DOO
 ☐ AOOR
 ☐ VOOR
 ☐ AAIR
 ☐ VVIR
 ☒ DOOR
 ☐ DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/>	Ventricular Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/>
Fixed AV Delay (ms)	<input type="text" value="70"/> <input type="text" value="70"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/>
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/>
Ventricular Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Activity Threshold (Low-High)	<input type="text" value="1"/> <input type="text" value="1"/>

Submit

User: zxc

Pacemaker Message

 Pacemaker Values Updated Successfully

OK

DDDR

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR **DDDR**

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="65"/> <input type="text" value="80"/> <input type="text" value="95"/> <input type="text" value="110"/> <input type="text" value="125"/> <input type="text" value="140"/> <input type="text" value="155"/> <input type="text" value="170"/>	Atrial Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="0.8"/> <input type="text" value="1.1"/> <input type="text" value="1.4"/> <input type="text" value="1.7"/> <input type="text" value="2.0"/> <input type="text" value="2.3"/> <input type="text" value="2.6"/> <input type="text" value="2.9"/> <input type="text" value="3.2"/>
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/> <input type="text" value="90"/> <input type="text" value="105"/> <input type="text" value="120"/> <input type="text" value="135"/> <input type="text" value="150"/> <input type="text" value="165"/>	Ventricular Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="0.8"/> <input type="text" value="1.1"/> <input type="text" value="1.4"/> <input type="text" value="1.7"/> <input type="text" value="2.0"/> <input type="text" value="2.3"/> <input type="text" value="2.6"/> <input type="text" value="2.9"/> <input type="text" value="3.2"/>
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/> <input type="text" value="60"/> <input type="text" value="70"/> <input type="text" value="80"/> <input type="text" value="90"/> <input type="text" value="100"/> <input type="text" value="110"/> <input type="text" value="120"/> <input type="text" value="130"/> <input type="text" value="140"/> <input type="text" value="150"/> <input type="text" value="160"/> <input type="text" value="170"/>	ARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/> <input type="text" value="180"/> <input type="text" value="210"/> <input type="text" value="240"/> <input type="text" value="270"/> <input type="text" value="300"/> <input type="text" value="330"/> <input type="text" value="360"/> <input type="text" value="390"/> <input type="text" value="420"/> <input type="text" value="450"/> <input type="text" value="480"/>
Fixed AV Delay (ms)	<input type="text" value="70"/> <input type="text" value="70"/> <input type="text" value="90"/> <input type="text" value="110"/> <input type="text" value="130"/> <input type="text" value="150"/> <input type="text" value="170"/> <input type="text" value="190"/> <input type="text" value="210"/> <input type="text" value="230"/> <input type="text" value="250"/> <input type="text" value="270"/> <input type="text" value="290"/>	VRP (ms)	<input type="text" value="150"/> <input type="text" value="150"/> <input type="text" value="180"/> <input type="text" value="210"/> <input type="text" value="240"/> <input type="text" value="270"/> <input type="text" value="300"/> <input type="text" value="330"/> <input type="text" value="360"/> <input type="text" value="390"/> <input type="text" value="420"/> <input type="text" value="450"/> <input type="text" value="480"/>
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	PVARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/> <input type="text" value="180"/> <input type="text" value="210"/> <input type="text" value="240"/> <input type="text" value="270"/> <input type="text" value="300"/> <input type="text" value="330"/> <input type="text" value="360"/> <input type="text" value="390"/> <input type="text" value="420"/> <input type="text" value="450"/> <input type="text" value="480"/>
Ventricular Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/> <input type="text" value="1.0"/> <input type="text" value="1.5"/> <input type="text" value="2.0"/> <input type="text" value="2.5"/> <input type="text" value="3.0"/> <input type="text" value="3.5"/> <input type="text" value="4.0"/> <input type="text" value="4.5"/> <input type="text" value="5.0"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="3"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="7"/> <input type="text" value="9"/> <input type="text" value="11"/> <input type="text" value="13"/> <input type="text" value="15"/> <input type="text" value="17"/> <input type="text" value="19"/> <input type="text" value="21"/> <input type="text" value="23"/> <input type="text" value="25"/> <input type="text" value="27"/> <input type="text" value="29"/>
Ventricular Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/> <input type="text" value="5"/> <input type="text" value="6"/> <input type="text" value="7"/> <input type="text" value="8"/> <input type="text" value="9"/> <input type="text" value="10"/> <input type="text" value="11"/> <input type="text" value="12"/> <input type="text" value="13"/> <input type="text" value="14"/> <input type="text" value="15"/> <input type="text" value="16"/>
Activity Threshold (Low-High)		<input type="text" value="1"/> <input type="text" value="1"/> <input type="text" value="2"/> <input type="text" value="3"/> <input type="text" value="4"/>	

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR **DDDR**

Lower Rate Limit (ppm)	50	Atrial Sense Threshold (V)	0.5
	50 65 80 95 110 125 140 155 170		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Upper Rate Limit (ppm)	75	Ventricular Sense Threshold (V)	0.5
	75 90 105 120 135 150 165		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Max Sensor Rate (ppm)	50	ARP (ms)	150
	50 60 70 80 90 100 110 120 130 140 150 160 170		150 180 210 240 270 300 330 360 390 420 450 480
Fixed AV Delay (ms)	70	VRP (ms)	150
	70 90 110 130 150 170 190 210 230 250 270 290		150 180 210 240 270 300 330 360 390 420 450 480
Atrial Amplitude (V)	0.5	PVARP (ms)	150
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		150 180 210 240 270 300 330 360 390 420 450 480
Ventricular Amplitude (V)	0.5	Reaction Time (s)	1
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Pulse Width (ms)	1	Recovery Time(s)	5
	1 2 3 4 5 6 7 8 9 10		5 7 9 11 13 15 17 19 21 23 25 27 29
Ventricular Pulse Width (ms)	1	Response Factor (Slow-Fast)	1
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)	1		
	1 2 3		

Submit

User: zxc Atrium Egram OK

Communication OK

i All parameters were transmitted properly

OK

Pacemaker Interface

☐ AOO
 ☐ VOO
 ☐ AAI
 ☐ VVI
 ☐ DOO
 ☐ AOOR
 ☐ VOOR
 ☐ AAIR
 ☐ VVIR
 ☐ DOOR
 ☒ DDDR

Lower Rate Limit (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	Atrial Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>
	50 65 80 95 110 125 140 155 170		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Upper Rate Limit (ppm)	<input type="text" value="75"/> <input type="text" value="75"/>	Ventricular Sense Threshold (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>
	75 90 105 120 135 150 165		0.5 0.8 1.1 1.4 1.7 2.0 2.3 2.6 2.9 3.2
Max Sensor Rate (ppm)	<input type="text" value="50"/> <input type="text" value="50"/>	ARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	50 60 70 80 90 100 110 120 130 140 150 160 170		150 180 210 240 270 300 330 360 390 420 450 480
Fixed AV Delay (ms)	<input type="text" value="70"/> <input type="text" value="70"/>	VRP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	70 90 110 130 150 170 190 210 230 250 270 290		150 180 210 240 270 300 330 360 390 420 450 480
Atrial Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	PVARP (ms)	<input type="text" value="150"/> <input type="text" value="150"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		150 180 210 240 270 300 330 360 390 420 450 480
Ventricular Amplitude (V)	<input type="text" value="0.5"/> <input type="text" value="0.5"/>	Reaction Time (s)	<input type="text" value="1"/> <input type="text" value="1"/>
	0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0		1 3 5 7 9 11 13 15 17 19 21 23 25 27 29
Atrial Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>	Recovery Time(s)	<input type="text" value="5"/> <input type="text" value="5"/>
	1 2 3 4 5 6 7 8 9 10		5 7 9 11 13 15 17 19 21 23 25 27 29
Ventricular Pulse Width (ms)	<input type="text" value="1"/> <input type="text" value="1"/>	Response Factor (Slow-Fast)	<input type="text" value="1"/> <input type="text" value="1"/>
	1 2 3 4 5 6 7 8 9 10		1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16
Activity Threshold (Low-High)		<input type="text" value="1"/> <input type="text" value="1"/>	
		1 2 3	
Submit			
User: zxc		Atrium Egram Ven	OK

Pacemaker Message

Pacemaker Values Updated Successfully

OK

4.4 Rate-Adaptive Modes:

The rate adaptive modes add some new parameters to the DCM interface. As a result, there are more possible test cases.

For any rate adaptive modes, the Max Sensor Rate should be greater than or equal to the Lower Rate Limit. In the case that it is not, a popup window that notifies the user and then prompts them to retry. In the below test case, we have such a situation, and when the values are submitted it performs as expected.

The screenshot displays the 'Pacemaker Interface' window. At the top, a mode selection bar includes options: AOO, VOO, AAI, VVI, DOO, **AOOR** (highlighted), VOOR, AAIR, VVIR, DOOR, and DDDR. Below this, several sliders and input fields are arranged in a grid:

- Lower Rate Limit (ppm):** A slider ranging from 50 to 170, with a current value of 80.
- Atrial Pulse Width (ms):** A slider ranging from 1 to 10, with a current value of 1.
- Upper Rate Limit (ppm):** A slider ranging from 75 to 165, with a current value of 140.
- Reaction Time (s):** A slider ranging from 1 to 29, with a current value of 1.
- Max Sensor Rate (ppm):** A slider ranging from 50 to 170, with a current value of 70.
- Recovery Time(s):** A slider ranging from 5 to 29, with a current value of 5.
- Atrial Amplitude (V):** A slider ranging from 0.5 to 5.0, with a current value of 0.5.
- Response Factor (Slow-Fast):** A slider ranging from 1 to 16, with a current value of 1.
- Activity Threshold (Low-High):** A slider ranging from 1 to 4, with a current value of 1.

At the bottom left, the text 'User: zxc' is displayed. In the center, there is a 'Submit' button. At the bottom right, there are three buttons: 'Atrium Egram', 'Ventricle Egram', and 'Dual Egram'.

Pacemaker Interface

☐ AOO
 ☐ VOO
 ☐ AAI
 ☐ VVI
 ☐ DOO
 ☐ AOOR
 ☐ VOOR
 ☐ AAIR
 ☐ VVIR
 ☐ DOOR
 ☐ DDDR

Lower Rate Limit (ppm)
 Atrial Pulse Width (ms)

Upper Rate Limit (ppm)
 Reaction Time (s)

Max Sensor Rate (ppm)
 Recovery Time(s)

Atrial Amplitude (V)
 Response Factor (Slow-Fast)

Activity Threshold (Low-High)

Submit

User: zxc

Pacemaker Message

i Max Sensor Rate must be greater than or equal to the Lower Limit.

OK

4.5 Improper Transmission/Other Issues

Some other issues can arise for each of the modes. For instance, if the Lower Rate Limit slider is above the Upper Rate Limit, there is a special case in the code to catch this. It sets the Lower Rate Limit equal to 50 ppm. Here is this case in action:

Pacemaker Interface

AOO VOO AAI VVI DOO AOOR VOOR AAIR VVIR DOOR DDDR

Lower Rate Limit (ppm) 114

Upper Rate Limit (ppm) 75

Ventricular Amplitude (V) 0.5

Ventricular Pulse Width (ms) 1

Submit

User: zxc

Atrium Egram Ventricle Egram Dual Egram

Pacemaker Message

Lower Rate Limit Fixed to 50ppm
Pacemaker Values Updated Successfully

OK